



KIT

Airkit Disposable Fold Flat Mask

Product Code

4026

Available Sizes

Mens, Ladies

Available Colours

White

Safety Standards

EN149

Packaging

20 / 400

Features

Hospitals, Medical Use, Veterinary

Applications

Hospitals, Medical Use, Veterinary



Currently out of stock due to Coronavirus outbreak*

The AIRKIT® particulate respirator 4026, EN149:2001 FFP2 provides protection against non-toxic solid and liquid aerosols in concentrations up to 12 x MAC*/OEL*/TLV* or 10 x APF in the UK. Disposable respiratory fold flat mask P2 with metal nose strip

Lightweight comfort mask with two piece comfortable headband, and adjustable noseband for fitting purposes. For protection against low - toxicity solids & liquids. P2 - For use with low to average toxicity particles e.g.

Cement works, hospitals, laboratories, mining, farming



KIT

Airkit Disposable Fold Flat Mask

Limitations of Use

For protection against low - toxicity solids & liquids. These respirators help reduce exposure to certain particles. Misuse may result in sickness or death. Before use, the wearer must read and understand the User Instructions provided as a part of the product packaging. Time use limitations may apply.

Storage Instructions

Vapours, including those present in paint spraying operations, asbestos, arsenic, cadmium, lead, methylenedianiline (MDA) or sand blasting. Concentrations of acid gases exceeding the OSHA PEL or applicable exposure limits, whichever is lower. Aerosol concentrations that exceed 10 times the OSHA PEL, or applicable exposure limits, whichever is lower. This respirator does not supply oxygen.

Maintenance Instructions

These respirators help reduce exposure to certain particles. Misuse may result in sickness or death. Before use, the wearer must read and understand the User Instructions provided as a part of the product packaging. Time use limitations may apply.