



# KIT

## Airkit Medifit Medical Mask FFP3

### Product Code

4333

### Available Sizes

Mens, Ladies L/XL S/M 4333S

### Available Colours

Amber, Blue

### Safety Standards

EN149

EN14683

### Packaging

10 / 200

### Features

Hospitals, Medical Use, Student, Veterinary

### Applications

Agriculture, Chemical Spraying, Construction, Councils, DIY, Factory, General Handling, Hospitals, Household Utility, Janitorial, Landscape Gardener, Local Authorities, Medical Use, Pharmaceutical, Student, Veterinary, Warehouse



\*\*\*Currently out of stock due to Coronavirus outbreak\*\*\*\*

This medical mask has a fluid resistant barrier. The Airkit particulate conical respirators FFP3 provide protection against non toxic solid and liquid in concentrations up to 50 x MAC/OEL/TDV or 20 APF in the UK..



# KIT

## **Airkit Medifit Medical Mask FFP3**

### **Limitations of Use**

For use with high toxicity & radioactive particles e.g. Pharmaceutical industry, toxic waste disposal, battery manufacture. These respirators help reduce exposure to certain particles. Misuse may result in sickness or death. Before use, the wearer must read and understand the User Instructions provided as a part of the product packaging. Time use limitations may apply. For proper use, see package instructions or consult you safety supervisor

### **Storage Instructions**

Vapours, including those present in paint spraying operations, asbestos, arsenic, cadmium, lead, methylenedianiline (MDA) or sand blasting. Concentrations of acid gases exceeding the OSHA PEL or applicable exposure limits, whichever is lower. Aerosol concentrations that exceed 10 times the OSHA PEL, or applicable exposure limits, whichever is lower. This respirator does not supply oxygen

### **Maintenance Instructions**

If the respirator becomes damaged, soiled or breathing becomes difficult, leave the contaminated area immediately and dispose of the respirator.